EASY WATER FAMILY

WATER EQUIVALENT PHANTOMS FOR QA IN RT







QUALITY ASSURANCE IN RADIOTHERAPY

Water-based and water-equivalent phantoms are qualified to measure dose distribution since the human body consits mainly of water. Water tank phantoms are used for acceptance, commissioning and annual QA tests. For frequently recurring tasks in quality assurance, water-equivalent RW3 material is an integral and time-efficient tool for every medical physicist.

Our EASY SLAB and EASY CUBE phantoms are made of RW3 and therefore perfectly suited for these requirements.

By its modular design and various extensions, the EASY CUBE phantom can also be used for socalled End-to-End QA tests to verify the whole treatment chain.

LAP'S EASY WATER FAMILY

- Highly homogenous material
- Easy to handle
- Various extensions and accessories
- Certified medical device

☑ EASY CUBE

MULTI-MODULAR PHANTOM FOR DOSIMETRY AND BEYOND

LAP IS YOUR IDEAL PARTNER

- Since 1984
- In-house hardware and software development
- Scientific collaborations
- Worldwide service network
- Certified in accordance to ISO 9001 and ISO 13485
- Made in Germany



EASY SLAB

SLAB PHANTOM FOR FREQUENT STANDARD DOSIMETRY





EASY SLAB

EASY SLAB supports frequently reccurent standard dosimetry in the field of radiation oncology.

The phantom is used in combination with irradiation detectors like ionization chambers. Corresponding adapter plates allow the detectors to be positioned as required for the QA task.

QA OF HIGH-ENERGY PHOTONS AND ELECTRONS ON LINACS

- Measurements with a homogeneous, standardized phantom with properties equivalent to water
- Measurements with irradiation detectors
- from all well-known vendors by use of optional adapter plates
- Measurement of dose distributions which have been computed for the phantom within the framework of radiation therapy planning
- Measurement of dose distributions which have been computed for the phantom in connection with intensitymodulated fields



EASY SLAB

EASY SLAB can be composed to a block phantom with a base of 300 mm \times 300 mm or 400 mm \times 400 mm and a maximum heigth of 300 mm. Detectors can be positioned in various measuring depths in steps of 1 mm using dedicated detector adapter plates.

For each of the two phantom sizes, a variety of different detector adapter plates are available.



■ is available as a cube-shaped phantom with an edge length of 300 mm

■ consists of 1 slab of 1 mm, 2 slabs of 2 mm, 1 slab of 5 mm and 29 slabs



EASY SLAB WILL BE SUPPLIED WITH A TRANSPORT AND STORAGE CASE.







THE UNIQUE MULTI-COLORED GRID EASES INITAL POSTIONING DURING SET UP AND ORIENTATION DURING QA TASKS AND THEREFORE REDUCES MISALIGNMENT RISKS.

EASY CUBE supports QA measurements in dosimetry and beyond in the field of radiation oncology. The phantom can be used in combination with irradiation detectors like ionization chambers, radiosensitive films and TLDs. The detectors are positioned inside the phantom with appropriate adapters which consist of the same material as the phantom.

QA OF HIGH-ENERGY PHOTONS AND ELECTRONS ON LINACS

THE FOLLOWING APPLICATIONS ARE POSSIBLE:

- Measurements with a homogeneous, standardized phantom with properties equivalent to water
- Measurements with irradiation detectors from all well-known vendors by use of optional adapters
- Measurement of dose distributions which have been computed for the phantom within the frame of radiation therapy planning
- Measurement of dose distributions which have been computed for the phantom in connection with intensity-modulated fields
- Measurement of dose distributions which have been computed for the phantom concerning stereotactic irradiation

MEASUREMENTS ON COMPUTER TOMOGRAPHS (CT)

THE FOLLOWING MEASUREMENTS ARE POSSIBLE:

- Cross-checking of Hounsfield values measured by a CT scanner
- Cross-checking of CT scans for geometric distortion
- Cross-checking of table coordinates displayed by a CT scanner
- Cross-checking of coordinate system of a CT scanner

FEATURES



BY MEANS OF VARIOUS EXTENSION MODULES, EASY CUBE CAN BE SHAPED INTO DIFFERENT PHANTOMS FOR DIFFERENT APPLICATIONS.

- Multi-modular design
- Different configurations possible
- Robust fixation of plates and inserts
- Reproducible measurement set-up



EASY CUBE WITH BODY MODULE **S**



EXTENSION MODULES

DIFFERENT SHAPES - DIFFERENT POSSIBILITIES



CYLINDER MODULE F OR S

The combination of EASY CUBE with the Cylinder Module **F** or Cylinder Module **S** forms a cylinder phantom with a diameter of 320 mm.

F: FILLABLE S: SOLID



HEAD / OVAL MODULE F

As with EASY CUBE, the Head / Oval Module ${\bf F}$ is hollow and can be filled with detectors and compensation elements.

This arrangement is especially suitable if radiation is to hit the phantom from different angles. It serves to prevent artefacts resulting from radiation hitting the edges of the phantom.

The Head / Oval Module **F** can be used to simulate a human head by fixating the two half cylinders to one another with assembly plates. The result is a cylindrical phantom, 180 mm in diameter and 180 mm in height.



BODY MODULE S

For simulation of the entire human torso (trunk) the Body Module **S** can be combined with the EASY CUBE. The Body Module **S** is comprised of two solid half cylinders and 16 body shaped distance plates with needles for film alignment (15 plates each 10 mm thick, numbered from 1 to 15, one plate 5 mm thick).



A MODULAR EQUIPMENT SOLUTION ACCESSORIES

INHOMOGENEOUS INSERTS

For calibration of Hounsfield values, various inhomogeneous inserts similar to fat, muscle, lung and bone fabrics are available. Artificial inhomogeneities like metal prostheses can be simulated as well.

The inhomogeneous inserts can be freely positioned inside EASY CUBE using the compensation elements.

DETECTOR ADAPTERS

- Adapters for different ionization chambers:
 - Adapter plates for radiochromic films: Radiochromic films can be placed inside EASY CUBE using compensation plates with a milled cavity. Closing of EASY CUBE, so it is lightproof, is not necessary.
 - Adapter plates for TLDs: or circular pits, so called "chips".

STEREOTACTIC LOCALIZER PLATES

Stereotactic localizer plates with integrated metal wires are available for the measurement of stereotactic coordinates. The stereotactic localizer plates are fastened with screws to the surface of the EASY CUBE.



COMPENSATION ELEMENTS

- Compensation plates of different thicknesses
- Compensation bars and spacers in different dimensions

LEVELING PLATES

Using the leveling plates, available in two different sizes, EASY CUBE and its extension modules can be fixed and precisely positioned within the application environment.



be delivered with a dedicated transport and storage case.







- One or more detector adapters corresponding to your configuration, e.g. for detectors already in your possession allow positioning at the desired measurement position.
- Different adapter plates are available for measurements with thermoluminescent dosimeters (TLD). The square adapter plates (160 mm edge length) are 10 mm thick. Depending on the TLD type, the adapter plates are equipped with bore holes, rods,





TRANSPORT AND STORAGE CASE

EASY CUBE together with the different modules of choice will



END-TO-END TEST

QA THROUGHOUT THE RADIATION TREATMENT CHAIN

With EASY CUBE we offer efficient QA from imaging to beam delivery which is cost-effective and fast.

EXAMPLES FOR QA WITH EASY CUBE

End-to-End tests need to be developed and coordinated closely to your individual clinical processes. Here are some options for your quality assurance througout your treatment chain.

QA OF CT PARAMETER

 \ni Starting with the alignment as preparation of the

⇒ Measurements of CT parameters like HU values

QA OF DATA TRANSFER

 \ominus As target volumes and OARs often are

 \ni A correct data transfer is mandatory

WHY END-TO-END TESTS?

End-to-end testing is a methodology used to test whether the treatment chain is performing as designed from start to finish. With the increasing complexity of external beam therapy, so-called End-to-End tests are intended to cover all steps from therapy planning to follow-up to fulfill the high demands on quality assurance. End-to-end tests are used to measure the overall accuracy of the radiation therapy chain, excluding patient specific factors. An End-to-End test is a prerequisite to the overall success of any IGRT to analyze potential errors accumulated by individual devices and processes in the treatment chain.

QA OF TREATMENT POSITIONING

QA OF TPS

⊖ In case of contouring and treatment planning are

⇒ An individual patient treatment plan can be

⇒ The EASY CUBE phantom can be aligned by use

⇒ Produced images can be checked regarding

⇒ Evaluated shifts can be checked to the realized

QA OF DOSE

⇒ Interesting values of dose distribution

⇒ Planned irradiation can be checked by





TECHNICAL DATA

Material	Water-equivalent white polystyrene "RW3" for high-energy photon and electron radiation				
		300 × 300 × 300 mm		1 slab of 1 mm thickness, 2 slabs of 2 mm, 1 slab of 5 mm, 29 slabs of 10 mm	
Energy range	Photons ⁶⁰ Co - 25 MV Electrons 4 - 23 MeV	400 × 400 × 300 mm		1 slab of 1 mm thickness, 2 slabs of 2 mm,	
Material	Polystyrene (C_8H_8) with admixture of			I slad of 3 mm, 29 slads of 10 mm	
composition	2.1% ± 0.2 % TiO ₂	EASY CUBE + Extension Modules			
Mass density	1.045 g /cm³	FASY CLIBE	180 x 180 x 180 mm		
(Z/A), value	0.536	LAST CODE			
EL 1 1		Cylinder Module		320 mm diameter, 180 mm height	
e/g)	3.386 × 10 ²³	Head / Oval		Head: 180 mm diameter, 180 mm height	
Electron concen- tration (e /cm³)	3.39 × 10 ²³	Iviodule	Oval: 300 mm lengin, 160 mm Width, 180 mm height		
		Body Module	360 mm length, 335 mm width, 180 mm height		





LAP GmbH

 Laser Applikationen

 Zeppelinstrasse 23

 21337 Lueneburg

 Germany

 Phone
 +49 4131 9511.95

 Fax
 +49 4131 9511.96

 Email
 info@lap-laser.com

www.lap-laser.com

LAP of America, LLC 161 Commerce Rd., Suite 3

Boynton Beach, FL 33426 USA Phone +1 561 416-9250 Fax +1 561 416-9263 Email america@lap-laser.com

LAP GmbH Laser Applikationen

 Представительство в Москве

 1, Казачий переулок 7

 119017 Москва

 Российская Федерация

 Тел.
 +7 495 7304043

 Факс
 +7 495 7304044

 Email
 info-russia.hc@lap-laser.com

LAP Laser Applications Asia Pacific Pte. Ltd.

750A Chai Chee Road #07-07 Viva Business Park Singapore 469001 Phone +65 6536 9990 Fax +65 6533 6697 Email info-asia.hc@lap-laser.com

LAP Laser Applications China Co. Ltd.

East Unit , 4F Building # 10 LujiaZui Software Park No. 61 Lane 91 EShan Road Shanghai 200127 China Phone +86 21 5047-8881 Fax +86 21 5047-8887 Email info-cn.hc@lap-laser.com

