

Optional accessories

ASTOR

Optional



Bracket for false ceilings

- Stable bracket for large distances between mounting plates and false ceiling or wall recess
- Length up to 1.500 mm




Free standing column

- For locations where wall mounts are not possible
- Length options: 900–1.400 mm or 1.100–1.600 mm



2nd swivel mount

- Bracket to enable 2nd tilting axis
- i.e. for mounting sagittal lasers at diagonally placed machines
- Swivel angle up to 45°



Mount for double floor

- Mount for rooms with framed or double floors
- Length up to 500 mm



Central power supply

- Enables central switching of up to 4 devices in one room
- Including one power supply, T-Connectors as needed, cables as needed
- Cable length from central switch to device: 15 m



WILKE Phantom

- For regular laser QA

LAP GmbH Laser Applikationen
 Zeppelinstr. 23
 21337 Lüneburg
 Germany

P +49 4131 95 11-95
 E info@lap-laser.com

 made
in
Germany