



Data Sheet

Levelling plate for THALES 3D SCANNER

Comfortable phantom positioning on the LINAC couch

Benefits

- Tilt corrections along transversal and longitudinal axis with high accuracy
- Efficient manual tilt alignment and visual levelling tool
- Engraved marker lines for alignment with room or machine lasers and incorporated positioning bolts for accurate phantom setup
- Compatible with most LINAC couches

Technical data

Dimensions (L × W × H)	1100 mm × 540 mm × 150 mm
Weight	18 kg
Material	High-resistance aluminium, 2 incorporated screw-on spirit levels
Compatibility	Compatible with both THALES versions (THALES 3D SCANNER and THALES 3D MR SCANNER)
Accuracy	Accuracy of the tilt correction better than $\pm 0.1^\circ$

Scope of delivery

- 1 levelling plate, assembled

LAP GmbH Laser Applikationen
Zeppelinstr. 23
21337 Lüneburg
Germany
P +49 4131 95 11-95
E info@lap-laser.com

